

# Centrifuge

**BUSINESS LINE:** Medical

**GROUP:** Centrifuge

**MEDICAL DEVICE**



## Accessories



**ROTOR\***

**MEDICAL DEVICE**



**BUCKETS**

**MEDICAL DEVICE**



**ADAPTORS**

**MEDICAL DEVICE**

## Consumables



**TUBE\***

**MEDICAL DEVICE**

### PURPOSE

Laboratory centrifuges are motor driven equipment, used to spin liquid samples at high speeds to separate the liquid samples in to their higher and lower densities. This works using the principle of sedimentation.

### GOOD TO KNOW

- The Rotors, Buckets and Adapters are selected for the specifics of the study and the tubes needed.
- Centrifuge Calibration is done yearly by sending a newly calibrated centrifuge to the site and returning the original for re-use.